



# Electrolytic Cell



## The Electrolytic Kit comes with

- 2 pieces of carbon rod
- 2 pieces of Pyrex (UK) test tube 16 x 150mm
- 26 cm wire with 2 crocodile clips at both ends

## To carry out Electrolysis Experiment

- Fill up the apparatus with solution (¾ full).
- Fill up test tubes (full), close it with your thumb
- Invert it under water, then into the carbon rod.
- Fit in the cover to hold the test tubes.
- Connect the crocodile clips to a 6.0-9.0 volt DC power supply

Overall Dimension of Electrolytic Kit (without test tube)	15cm (H) x 4cm (L) x 7cm (W) x 10.5cm (D)
Quantity per carton	10 pieces
Material	PS Plastic
Weight per piece	160gm
Net weight per carton	2.01kg
Overall Dimension per carton	380mm(L) x 170mm(W) x 197mm(H)
Country of Origin	<b>Malaysia</b>

**Important :** Always hold the body of the apparatus to avoid wobbling!

D - Depth    W - Width    L - Length    H - Height

We reserve the right to revise the price and colour of this product at any time in future without prior notice.

